

ENGLISH

Instructions For Use

WavelinQ™ EndoAVF System (WavelinQ™ System, WavelinQ™ or REF WQ4305) Components

- WavelinQ™ Arterial Catheter
- WavelinQ™ Venous Catheter

WavelinQ™ EndoAVF System Rated Accessory Voltage

- WavelinQ™ Venous Catheter is rated for 700Vpeak

Commercially Available Devices

- TVA or BD ESU-1 Electrosurgical Unit (ESU)
- Electrosurgical Pencil:
 - Must have a 3-prong universal connector that interfaces with the ESU-1 Electrosurgical Generator
 - Must have a minimum Rated Accessory Voltage of 700V and a Yellow activation button.
 - Ground Pad
 - Must have a universal connector that interfaces with the ESU-1 Generator
 - Must be rated to disperse a minimum of 60 W for 1 Second.
 - TZ Medical Arm Board and Fixation Straps: CZ-400-TVA (Arm Restraint) and TVA-MC-2 (Fixation Straps)

System Description

The WavelinQ™ EndoAVF System consists of two 4 Fr single-use, disposable, magnetic, hydrophilic coated catheters; a venous catheter and an arterial catheter. It is used with the ESU-1 Electrosurgical Unit and an Electrosurgical Pencil. The WavelinQ™ venous catheter is a flexible 4 Fr magnetic catheter that contains a radiofrequency (RF) electrode, a hemostasis valve crosser for interfacing with the introducer sheath on the distal end, and has a cable and plug for the delivery of radiofrequency energy. The WavelinQ™ venous catheter connects via the Electrosurgical Pencil to the ESU for delivery of radiofrequency energy with standard grounding pads. The WavelinQ™ arterial catheter is a flexible 4 Fr magnetic catheter that contains a backstop for receiving the electrode. When placed in proximity, the magnets contained in each catheter attract to each other, while aligning the electrode with the backstop. Rotational indicators are present in each catheter and used to accurately position the catheters. Radiofrequency energy can then be delivered through the electrode for cutting and/or coagulating tissue.




Figure 1. WavelinQ™ Venous and Arterial Catheter Cross-Section



Figure 2. WavelinQ™ Venous Catheter with Hemostasis Valve Crosser in place




Figure 3. WavelinQ™ Packaging and Catheter Identification

Table 1. Specifications	WavelinQ™ Venous Catheter	WavelinQ™ Arterial Catheter
Maximum guidewire size	0.014 in (0.356 mm)	0.014 in (0.356 mm)
Catheter working length, cm/in	42 cm (16.4 in)	50 cm (19.7 in)
Maximum diagonal of square section	1.55 mm (0.061 in)	1.55 mm (0.061 in)
Minimum sheath introducer inner diameter	1.78 mm (0.070 in)	1.78 mm (0.070 in)
Recommended sheath introducer size	5 Fr	5 Fr
Rated Accessory Voltage	700 Vp	-
Electrosurgical Generator Settings	Cut T, 60 W, 0.7 Sec	-

Indications

The WavelinQ™ EndoAVF System is intended for the cutting and coagulation of blood vessel tissue in the peripheral vasculature for the creation of an arteriovenous fistula used for hemodialysis.

Contraindications

- Known central venous stenosis or upper extremity venous occlusion on the same side as the planned AVF creation.
- Known allergy or reaction to any drugs/fluids used in this procedure.
- Known adverse effects to moderate sedation and/or anesthesia.
- Distance between target artery and vein > 1.5 mm.
- Target vessels < 2 mm in diameter.

Warnings, Cautions, and Precautions

Warnings

1. **The WavelinQ™ EndoAVF System is only to be used with the approved commercially available devices specified above. Do not attempt to substitute non-approved devices or use any component of this system with any other medical device system.**
2. **The WavelinQ™ catheters are single use devices. DO NOT re-sterilize or re-use either catheter. Potential hazards of reuse include infection, device mechanical failure, or electrical failure, potentially resulting in serious injury or death.**
3. **Use caution when performing electrosurgery in the presence of pacemakers.**
4. **Improper use could damage insulation that may result in injury to the patient or operating room personnel.**
5. **Do not plug device into the electrosurgical pencil with ESU on.**
6. **Keep active accessories away from patient when not in use.**
7. **Do not permit cable to be parallel to and/or in close proximity to leads of other devices.**
8. **Do not wrap cable around handles of metallic objects such as hemostats.**
9. **Consult the ESU User's Guide on its proper operation prior to use.**
10. **Do not use closure devices not indicated to close the artery used for access.**

Cautions

1. Only physicians trained and experienced in endovascular techniques should use the device.
2. Adhere to universal precautions when utilizing the device.
3. Do not kink, pinch, cut, bend, twist, or pull excessively or with excessive force on any portion of the devices. Damage to the catheter body may cause the device to become inoperable.
4. Avoid sharp bends. This may cause the device to become inoperable.
5. Do not pinch or grasp the catheter with excessive force or with other instruments. This may cause the device to become inoperable.
6. Do not bend the rigid portion of the catheter near the electrode or backstop.
7. Do not touch or handle the active electrode. Electrode dislodgement may occur.
8. Always use the hemostasis valve crosser to assist insertion of the venous catheter through the introducer sheath. Insertion into introducer sheath without hemostasis valve crosser may damage electrode.
9. Do not attempt to remove the hemostasis valve crosser located on the venous device. Device damage or fracture may occur.

Precautions

1. Care should be taken to avoid the presence of fluid on the ESU.
2. Care should be taken during handling of the arterial and venous catheters in patients with implantable cardiac defibrillators or cardiac pacemakers to keep the distal 3 inches of the catheters at least 2 inches from the implanted defibrillator or pacemaker.
3. Care should be taken to avoid attempting fistula creation in a heavily calcified location of a vessel as fistula may not be adequately formed.
4. The safety and performance of this device has not been established for pediatric patients.
5. If the device does not perform properly during the creation of the endovascular fistula it is possible that a fistula will not be created or there may be some vessel injury.
6. Keep magnetic ends of catheters away from other metallic objects which may become attracted and collide with devices.

Electrical Safety

The WavelinQ™ EndoAVF System complies with the electrical safety standards 60601-1, 60601-1-2 and 60601-2-2.

Potential Adverse Events

The known potential risks related to the WavelinQ™ EndoAVF System and procedure, a standard AVF, and endovascular procedures may include, but are not limited to: aborted or longer procedure; additional procedures; bleeding, hematoma or hemorrhage; bruising; burns; death; electrocution; embolism; failure to mature; fever; increased risk of congestive heart failure; infection; numbness, tingling, and/or coolness; occlusion/stenosis; problem due to sedation or anesthesia; pseudoaneurysm; aneurysm; sepsis; steal syndrome or ischemia; swelling, irritation, or pain; thrombosis; toxic or allergic reaction; venous hypertension (arm swelling); vessel, nerve, or AVF damage or rupture; wound problem.

How Supplied

Packaging

The WavelinQ™ EndoAVF System is supplied STERILE and intended for single use only. Catheters are considered sterile only if the pouch is undamaged and unopened. The outer surfaces of the carton and pouch are NON-STERILE and must not be placed in the sterile field. Open the pouch using aseptic techniques so that its inner contents are delivered to the sterile field. The WavelinQ™ catheters are sterilized with ethylene oxide gas.

Storage

Store the WavelinQ™ EndoAVF System in a cool, dark, dry place.

Additional Equipment

Note: While exhaustive, this equipment list is not meant to cover all possible scenarios.

Additional Equipment

- Ultrasound machine and ultrasound probe
- Two (2) 5 Fr introducer sheaths
- Two (2) 0.014" (0.356 mm) guidewires
- 4 Fr guide catheter (as needed)
- Tourniquet or blood pressure cuff
- Anticoagulant medication, as needed
- Vasodilator medication, as needed

Electrode Activation

- ESU-1 Electrosurgical unit
- Electrode pencil
- Ground pad
- Arm restraint and fixation straps

Instructions For Use

Pre-Procedural Preparations

1. The procedure should be performed in an angiography room and carried out under X-ray control.
2. Patient preparation and sterile precautions should be the same as for any percutaneous transcatheter procedure. The medication is decided by the physician, including anesthesia and precautions to reduce pain, clotting, and vasospasm during the procedure according to latest scientific guidelines and with respect to the individual patient.
3. Place ESU on a flat secure surface located near operative field making sure that the ground pad and electrosurgical pencil cabling have sufficient length to be connected during subsequent procedural steps.
4. Turn on ESU. Ensure the **Cut T** mode is illuminated, the power setting LED display reads **60 W**, and the maximum activation time of **0.7 SEC** is set in the time LED display.
5. Place ground pad on patient following standard guidelines for electrosurgical patient grounding and insert ground pad plug in to ESU. Confirm indicator light changes from red to green, ensuring appropriate patient contact.
6. Turn off ESU until ready for energy delivery in order to prevent inadvertent activation.

Vascular Access

7. Administer anesthesia or conscious sedation per hospital protocol.
8. Secure the patient's procedure arm in a restraint to prevent arm movement.
9. Use tourniquet or blood pressure cuff to facilitate vessel access as required per standard protocol.
10. Gain percutaneous access to target vein with a puncture needle.
11. Introduce a guidewire into the vein through the needle and advance an adequate length to facilitate introducer sheath insertion. A microintroducer sheath may be used to first confirm intraluminal position of guidewire.
12. Insert a 5 Fr sheath into the vein over the guidewire.
13. Inject contrast media and perform a venogram to assess appropriateness of vessel anatomy for procedure. Limit the dose of contrast depending on the patient's residual renal function. Alternate contrast methods may be used at the discretion of physician.
14. Gain access to target artery with a puncture needle and introduce a guidewire into the artery.
15. Insert a 5 Fr sheath into the artery over the guidewire.
16. Using physician discretion, administer anticoagulant and vasodilator intra-arterially or intravenously.
17. Insert a 0.014" guidewire into the artery and deliver to the target AV creation site.
18. Insert a 0.014" guidewire into the vein and deliver to the target AV creation site.
19. Orient fluoroscope perpendicular to the target artery and vein using guidewires, ultrasound, or contrast as guidance.
20. Remove tourniquet or blood pressure cuff (if applied).

Preparation of the Catheters

21. Carefully inspect the WavelinQ™ EndoAVF System pouch for any evidence of damage to the sterile barrier. If there is evidence of damage, do not use the WavelinQ™ EndoAVF System.
22. After inspection of the pouch, carefully peel open the pouch and transfer the sterile WavelinQ™ EndoAVF System to sterile field using standard transfer precautions.

EndoAVF Creation Procedure

23. Remove the arterial catheter from the packing card and inspect for damage (see Figures 1-3). Evaluate the distal end of the catheter. If it is suspected that the sterility or performance of the catheter has been compromised, the catheters should not be used.
24. Advance the arterial catheter over the 0.014" wire and insert through the arterial sheath, taking care not to kink the distal magnet arrays. Under fluoroscopic guidance, advance the catheter to the target AVF location.
25. Rotate the arterial catheter until the illumination of the rotational indicators is maximized and the concave surface of the backstop is pointed at the venous wire.
26. Remove the venous catheter from the packing card:
 - a) Removing the plug and cable bundle from the card tabs
 - b) Uncouple the venous handle and then slide the catheter out of its containment tube. See Figures 1-3.
27. Inspect the venous catheter for damage. If it is suspected that the sterility or performance of the catheter has been compromised, the catheter should not be used.
28. Catheters have a hydrophilic coating and at the physician's discretion may be hydrated prior to insertion by wiping with saline in the sterile field.
29. Advance the venous catheter over the 0.014" wire until the yellow hemostasis valve crosser encounters the hemostasis valve of the introducer sheath. Grasp and insert the yellow hemostasis valve crosser through the hemostasis valve until it stops in the sheath hub. Grasp the proximal end of the yellow valve crosser and advance the catheter simultaneously through the yellow crosser and the sheath. Fluoroscopically inspect the electrode after venous catheter insertion to confirm proper electrode form. If the electrode appears deformed, gently remove venous catheter and inspect. If upon direct visual inspection the electrode appears damaged, replace the venous catheter. If venous catheter is removed and requires reinsertion, reposition the yellow valve crosser over the electrode prior to advancing through the hemostasis valve.




Figure 4. Venous Catheter & Valve Crossover, Pre-insertion




Figure 5. Venous Catheter & Valve Crossover, Post-insertion

30. Under fluoroscopic guidance, advance the venous catheter towards the target location until the distal magnets of the venous catheter begin to engage the first proximal magnets of the arterial catheter. In this location, pause to rotate the venous catheter until the illumination of the rotational indicators are maximized and the arc of the electrode is pointed at the arterial catheter. Adjust the arterial catheter as needed to confirm that catheters are aligned and prepared for venous catheter advancement.
31. Advance the venous catheter until the arc of the electrode is congruent with concave surface of the arterial backstop. Electrode should appear compressed. Confirm that the rotational indicators appear aligned and rotationally similar.
32. Rotate the fluoroscope to visualize the maximum tissue thickness distance between arterial and venous catheters. Confirm that the tissue thickness adjacent to the electrode housing is no greater than the width of the magnet array which is 1mm. If tissue thickness appears greater, adjust catheter position to a thinner tissue segment.
33. With catheters in confirmed activation position, remove cable tie, remove preassembled insert, and then connect venous catheter plug to electro-surgical pencil. Fully insert plug until there is no metallic surface exposed.
34. Pass electro-surgical pencil hand switch connector out of the sterile field and connect it to ESU's Monopolar 1 receiver.
35. Retract or remove both 0.014" guidewires from the catheter activation zone. No guidewires should be present between the proximal and distal magnet zones during activation.
36. Ensure tourniquet has been removed (if applied).
37. Do not allow catheters to move in order to minimize chance of misalignment.
38. Using fluoroscopy, verify final catheter and electrode position before energy delivery.



Figure 6. Activation Position

39. Hold patient's procedure arm at the wrist with firm pressure to minimize arm flexion and rotation during energy delivery.
40. Turn on ESU and again ensure the **Cut T** mode is illuminated, the power setting LED display reads **60 W**, and the maximum activation time of **0.7 SEC** is set in the time LED display.
41. While imaging with fluoroscopy, deliver RF energy by firmly pressing and holding the **yellow** cut switch on the electro-surgical pencil until the audible ESU activation tone stops. The electrode should visibly advance and touch the arterial backstop. If electrode does not contact backstop, an additional activation may be administered under the condition that the catheters have not been moved from their original position. Do not activate the device more than 3 times.




Figure 7. Electrode Advanced to Backstop

42. Remove the venous catheter. Remove the arterial catheter.
43. Perform arteriogram via the arterial sheath to confirm AVF creation. Alternate contrast methods may be used at the discretion of physician.
44. Embolize a brachial vein if patient anatomy allows.
45. Remove the arterial sheath and venous sheath and complete the puncture site closure to achieve hemostasis.

Device Disposal

WaveInQ™ EndoAVF System components that have been in contact with body fluids are a potential biohazard. Handle and dispose of the catheters and components using acceptable medical practice and all applicable local, provincial, and federal laws and regulations.

WavelinQ™ EndoAVF System has been previously referred to as the everlinQ™ endoAVF System.

Summary of Supplementary Clinical Information (WavelinQ™ EndoAVF System)

Objective

The WavelinQ™ EndoAVF System was supported by a global analysis ("4Fr Global Analysis") that was designed to aggregate and analyze safety and performance data from the EASE, EASE-2, and EU Study studies of subjects. The three studies enrolled subjects with chronic kidney disease who need hemodialysis and were candidates for percutaneous AVF creation with the WavelinQ™ EndoAVF System.

Design

Safety and performance data were collected on all study subjects during the follow-up period and the results reported here are the 6-month results of all available data. The pooled population evaluated for the aggregate study was 91 subjects.

Safety Endpoints

The safety endpoints are defined as the proportion of subjects with the following:

1. Significant Events
2. Serious Adverse Events (SAEs)
3. Device-related SAEs
4. Procedure-related SAEs
5. Closure device-related SAEs
6. Coil-related SAEs

Effectiveness Endpoints

The effectiveness endpoints are defined as the proportion of subjects with the following:

1. Procedure Success: Successful endoAVF creation is confirmed via intraprocedural angiography/fistulogram or duplex ultrasound verification performed post-procedure. This definition corresponds to the term "Technical Success" used by others.
2. Time to cannulation: The interval of time from the index procedure to the first successful 2-needle cannulation of the endoAVF.
3. Cannulation success: A successful cannulation of the endoAVF with 2-needles for dialysis. A subject may be called a 'cannulation success' with the first successful 2-needle cannulation of the endoAVF.
4. Primary patency: The interval from the time of access placement until any intervention designed to maintain or reestablish patency, access thrombosis, access abandonment, or the time of measurement of patency (SVS Reporting Standards definition).
5. Modified primary patency: Identical to Primary Patency except that loss was also triggered by reinterventions not directly related to the access circuit; namely coiling or vessel ligation of venous outflow tributaries to encourage flow into the superficial, more easily accessible veins of the upper arm.
6. Assisted primary patency: The interval from access placement to thrombosis or abandonment; not triggered by access circuit interventions performed in the absence of occlusion.
7. Secondary patency: The interval from the time of access placement until access abandonment, lost to thrombosis, or the time of patency measurement including intervening manipulations (surgical or endovascular interventions) designed to re establish functionality in thrombosed access (SVS Reporting Standards definition).
8. Functional patency: The interval of time from the first 2-needle dialysis utilizing the access until access abandonment (SVS Reporting Standards definition).
9. Functional cannulation: Successful 2-needle access of the endoAVF access circuit with performance of more than 2/3rds of dialysis sessions of at least 120 minutes in duration over a continuous 28-day period. This measure was defined to more aptly measure whether an endoAVF resulted in a working access site for a subject, and is more pertinent to a subject, as opposed to Successful Cannulation, which is limited to completing a single successful 2-needle cannulation.

Safety Results

In total, 22.0% (20/91) of the 4 Fr endoAVF subjects experienced a significant event, defined as device or procedure-related adverse events that either could be limb-threatening if not promptly identified or treated, or required additional therapy to reestablish patency of the endoAVF access circuit, irrespective of whether they met the criteria for an SAE. These significant events included stenosis, occlusion, thrombosis, and pseudoaneurysm of the access circuit and/or endoAVF, as well as one subject with abandonment of the endoAVF after a cannulation induced brachial artery injury.

SAEs were reported in 26.4% (24/91) of the Pooled Population. Of these events, 18/24 (75.0%) were unrelated to the device and unrelated to the procedure. There were 3 device-related SAEs reported in the studies (3/91, 3.3%) and included 1 thrombosis of the endoAVF, 1 stenosis of the endoAVF, and 1 access circuit false aneurysm. There were 5 procedure related SAEs reported in 5.5% (5/91) of the Pooled Population. None of the SAEs were related to the method of arterial access for the procedure. There were no reports of Unanticipated Adverse Device Effects (UADEs). There were no closure device-related SAE's or coil related SAE's reported.

Effectiveness Results

Procedural success, defined as the successful creation of an endoAVF with blood flow confirmed intraoperatively by fistulography or postoperative duplex ultrasonography, was achieved in 96.7% (88/91) of the Pooled Population. Using Kaplan-Meier (K-M) point estimates at month 6, the cannulation success rate was 84.9% ± 5.3% for subjects in the dialysis subset. The mean time to cannulation was 2.0 ± 1.6 months. Functional cannulation rate was achieved in 91.7% ± 7.8% of dialysis subjects.

Primary, Assisted Primary, and Secondary patency were determined using the K-M point estimates at month 6. Primary Patency was achieved in 72.4% ± 5.2%. at 6 months. Corresponding rates for Assisted Primary and Secondary Patency were each 77.3% ± 5.0%. Modified Primary Patency was achieved in 69.7% (±5.4%) and Functional Patency was achieved in 100% at 6 months in the Pooled Population.

Conclusion

The aggregate data presented confirm the safety and performance of the WavelinQ™ EndoAVF System for the creation of AVF for hemodialysis access.

REF

Catalogue Number
 Numéro de catalogue
 Katalognummer
 Numero di catalogo
 Número de catálogo
 Catalogusnummer
 Número do catálogo
 Αριθμός καταλόγου
 Katalognummer
 Artikelnummer
 Luettelonumero
 Katalognummer
 Numer katalogowy
 Katalógusszám
 Katalogové číslo
 Katalog Numarasi
 目錄編號
 카탈로그 No.
 Номер по каталогу

LOT

Lot Number
 Numéro de lot
 Los Nummer
 Numero di lotto
 Número de lote
 Lotnummer
 Número do lote
 Αριθμός ποτίδας
 Lotnummer
 Lot-nummer
 Eränumero
 Lotnummer
 Numer serii
 Tételszám
 Číslo šarže
 Parça Numarasi
 批號
 로트 번호(Lot No.)
 Номер партии



Use By Date
 Date limite d'utilisation
 Verwendbar bis
 Utilizzare entro
 Usar antes de
 Te gebruiken vóór
 Prazo de validade
 Ημερομηνία λήξης
 Anvendes før
 Utgångsdag
 Käytettävä ennen
 Brukes innen
 Termin ważności
 Felhasználható
 Datum použitelnosti
 Son Kullanım Tarihi
 有效期限
 유효기한
 Срок годности



Do Not Re-use
 Ne pas réutiliser
 Nicht wiederverwenden
 Non riutilizzare
 No reutilizar
 Niet opnieuw gebruiken
 Não reutilizar
 Μην επαναχρησιμοποιείτε
 Ikke til fjerdegangbrug
 Får ej återanvändas
 Ei saa käyttää uudestaan
 Ikke til gjenbruk
 Nie używać ponownie
 Újrahasználni tilos
 K jednorázovému použití
 Tekrar Kullanmayiniz
 請勿重複使用
 재사용하지 마십시오
 Не подлежит повторному использованию



Do Not Resterilize
 Ne pas restériliser
 Nicht resterilisieren
 Non risterilizzare
 No reesterilizar
 Niet opnieuw steriliseren
 Não reesterilizar
 Μην επαναποστειρωσετε
 Má ikke resteriliseres
 Får ej omsteriliseras
 Ei saa steriloida uudestaan
 Má ikke resteriliseres
 Nie sterylizować ponownie
 Újrasterilizálni tilos
 Neprovádějte resterilizaci
 Tekrar Sterilize Etmeyiniz
 請勿重複消毒
 재멸균하지 마십시오
 Повторная стерилизация запрещается



Contents
 Contenu
 Inhalt
 Contenuto
 Contenido
 Inhoud
 Conteúdo
 Περιεχόμενα
 Inndhold
 Innehåll
 Sisältö
 Innhold
 Zawartość
 Tartalom
 Obsah
 İçerikler
 炭物
 내용물
 Комплектация



Consult Instructions for Use
 Consulter le mode d'emploi
 Gebrauchsanweisung beachten
 Leggere le istruzioni per l'uso
 Consulte las Instrucciones de uso
 Raadpleeg de gebruiksaanwijzing
 Consultar as instruções de utilização
 Συμβουλευτείτε τις οδηγίες χρήσης
 Se bruksanvisningen
 Se bruksanvisning
 Lue käyttöohjeet
 Se bruksanvisning
 Zapoznać się z instrukcją obsługi
 Lásd a használati útmutatót
 Řiďte se návodem k použití
 Kullanna Talimatlarına Başvurun
 請參考使用說明
 사용 지침 참조
 Обратитесь к инструкциям по применению



Keep Away From Sunlight
 Conserver à l'abri de la lumière du soleil
 Von Sonneneinstrahlung fern halten
 Tenere al riparo dalla luce solare
 Manterer alejado de la luz solar
 Uit de buurt van zonlicht houden
 Manter afastado da luz solar
 Φυλάσσετε το προϊόν μακριά από το ηλιακό φως
 Beskyttes mod sollys
 Får ej utsättas för solljus
 Suojattava auringon valolta
 Beskyttes mot sollys
 Chronić przed promieniowaniem słonecznym
 Napfénytől védve tartandó
 Nevystavujte púsobení slnečného záření
 Güneş Işığından Uzak Tutunuz
 避免日光照射
 직사광선에 노출시키지 마십시오
 Оберегать от воздействия солнечных лучей



Keep Dry
 Conserver à l'abri de l'humidité
 Trocken halten
 Mantere asciutto
 Manter seco
 Droog bewaren
 Manter seco
 Διατηρείτε το προϊόν στεγνό
 Opbevares tørt
 Förvaras torr
 Suojattava kosteudelta
 Oppbevares tørt
 Chronić przed wilgocią
 Szárazon tartandó
 Uchovávejte v suchu
 Kuru Tutunuz
 保持乾燥
 건조한 상태로 보관
 Хранить в сухом месте



Not Made with Natural Rubber Latex
 Fabriqué sans latex naturel
 Nicht aus Naturgummitex hergestellt
 Non prodotto con lattice di gomma naturale Este producto no se fabrica con látex de cacho natural
 Niet gemaakt met latex van natuurlijk rubber
 Não contém látex de borracha natural
 Δεν κατασκευάζεται από φυσικό ελαστικό λάτεξ
 Ikke fremstillet med naturlig gummitæx
 Ej tilbered med naturgummitæx
 valmistuksessa ei ole käytetty luonnonkummitæksia
 Dette produktet er ikke produceret med naturlig gummitæks
 Produkt nie jest wytworzony z lateksu kauczuku naturalnego
 Nem tartalmaz természetes gummitæxet
 Není vyroben z přírodního pryžového latexu
 Doğal Kaucuk Latekten Üretilmemiştir
 製造未採用天然膠乳
 천연 고무 라텍스를 사용해 제조되지 않음
 Не содержит натуральный каучуковый латекс

**Do Not Use if Package is Damaged**

Ne pas utiliser si l'emballage est endommagé
 Nicht verwenden, wenn die Verpackung beschädigt ist
 Non utilizzare se la confezione è danneggiata
 No usar si el envase está dañado
 Niet gebruiken wanneer de verpakking beschadigd is
 Não utilizar se a embalagem estiver danificada
 Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία του έχει υποστεί ζημιά
 Má ikke bruges, hvis emballagen er beskadiget
 Använd inte produkten om förpackningen skadats
 Ei saa käyttää, jos pakkaus on vaurioitunut
 Skal ikke brukes hvis pakningen er skadet
 Nie stosować, jeśli opakowanie jest uszkodzone
 Ne használja fel, ha a csomagolás sérült
 Nepoužívejte, pokud je obal poškozen
 Paket Hasarlıysa Kullanmayınız
 如果包裝已受損，請勿使用
 포장에 손상이 있다면 사용하지 마십시오
 Не использовать, если упаковка повреждена

**Maximum Guidewire**

Guide maximal
 Maximalgröße Führungsdraht
 Dimensione massima del filo guida
 Máximo de guía
 Maximale geleidedraad
 Fio-guia máximo
 Μέγιστο όδηγό σύρμα
 Maks. guidewire
 Max ledare
 Ohjainvaijeri enintään
 Maksimal ledevaier
 Maksymalny rozmiar prowadnika
 Maximális vezetődírt
 Maximální vodič drát
 Maksimum Kılavuz Tel
 最大導線
 최대 가이드와이어
 Максимальный диаметр проводника

**Sterilized By Using Ethylene Oxide**

Stérilisé à l'oxyde d'éthylène
 Sterilisiert mit Ethylenoxid
 Sterilizzato mediante ossido di etilene
 Esterilizado mediante óxido de etileno
 Gesteriliseerd met behulp van ethyleenoxide
 Esterilizado por óxido de etileno
 Αποστειρωμένο με αιθυλενοξειδίο
 Steriliseret ved etylenoxid
 Steriliserad med etylenoxid
 Steriloitu etyleenoksidilla
 Steriliseret med etylenoksid
 Produkt sterylizowany tlenkiem etylenu
 Etilén-oxidál sterilizálva
 Etilén-oxidál etylenoxidem
 Etilen Oksit ile Sterilize Edilmıştır
 使用環氧乙烷消毒
 산화에틸렌 멸균
 Стерилизовано этиленоксидом

**Minimum Introducer**

Introduceur minimum
 Mindestgröße Einführhilfe
 Introduttore minimo
 Introdutor mínimo
 Minimale inbrenger
 Introdutor mínimo
 Ελάχιστος εισαγωγέας
 Mindste introducer
 Minsta införingsanordning
 Pienin sisäänviejä
 Minste innfører
 Minimalny rozmiar introduktora
 Minimálisan szűkséges bevezető
 Minimální zavadeč
 Minimum Introduser
 最小導引器
 최소 유도관
 Минимальный калибр интродьюсера

**Non-Pyrogenic**

Apyrogène
 Pyrogenfrei
 Apirogeno
 Apirógeno
 Niet-pyrogeen
 Apyrogénico
 Μη πυρογόνο
 Pyrogenfri
 Pyrogenfri
 Pyrogeeniton
 Pyrogenfri
 Apyrogenny
 Pirogénmentes
 Apyrogéní
 Pirojenik degildir
 無熱原
 비발열성
 Апиrogenно

**External Diameter**

Diamètre externe
 Außendurchmesser
 Diametro esterno
 Diámetro externo
 Buitendiameter
 Diámetro externo
 Ὑδενδிக diameter
 Yttre diameter
 Ulkoläpimitta
 Utvendig diameter
 Średnica zewnętrzna
 Külső átmérő
 Vnější průměr
 Dış Çap
 外徑
 외경
 Внешний диаметр

**Manufacturer**

Fabricant
 Hersteller
 Produttore
 Fabricante
 Fabricante
 Κατασκευαστής
 Producent
 Tilverkare
 Valmistaja
 Produzent
 Producent
 Gyártó
 Výrobce
 Üretici
 製造商
 제조사
 Производитель

**Type B Applied Part**

Pièce de type B
 Anwendungsteil vom Typ B
 Componente applicato di tipo B
 Parte aplicada de tipo B
 Type B toegepast onderdeel
 Peça aplicada de tipo B
 Εξάρτημα εφαρμογής τύπου ΒF
 Type B anvendt del
 Typ B-tillämpad del
 Typin B liityntäosa
 Anvendt del av type B
 Część aplikacyjna typu B
 B típusú alkalmazott alkatrész
 Použitá část typu B
 B Tipi Uygulamalı Parça
 B 型觸身部件
 유형 B 적용 부품
 Контактующая часть типа B

**Date of Manufacture**

Date de fabrication
 Herstellungsdatum
 Data di fabbricazione
 Fecha de fabricación
 Productiedatum
 Data de fabrico
 Ημερομηνία κατασκευής
 Fremstillingsdato
 Tilverningsdatum
 Valmistuspäivämäärä
 Produksjonsdato
 Data produkcji
 Gyártás dátuma
 Datum výroby
 Üretim Tarihi
 製造日期
 제조 연월일
 Дата производства



Copyright © 2020 Becton, Dickinson and Company. All rights reserved

Copyright © 2020 Becton, Dickinson and Company. Tous droits réservés.
 Copyright © 2020 Becton, Dickinson and Company. Alle Rechte vorbehalten.
 Copyright © 2020 Becton, Dickinson and Company. Tutti i diritti riservati.
 Copyright © 2020 Becton, Dickinson and Company. Todos los derechos reservados.
 Copyright © 2020 Becton, Dickinson and Company. Alle rechten voorbehouden.
 Copyright © 2020 Becton, Dickinson and Company. Todos os direitos reservados.
 Copyright © 2020 Becton, Dickinson and Company. Με την επιφύλαξη παντός δικαιώματος.
 Copyright © 2020 Becton, Dickinson and Company. Alle rettigheder forbeholdes.
 Copyright © 2020 Becton, Dickinson och Company. Eftertryck förbjudes.
 Tekijänoikeudet © 2020 Becton, Dickinson and Company. Kaikki oikeudet pidätetään.
 Copyright © 2020 Becton, Dickinson and Company. Alle rettigheter forbeholdt.
 Copyright © 2020 Becton, Dickinson and Company. Wszelkie prawa zastrzeżone.
 Copyright © 2020 Becton, Dickinson and Company. Minden jog fenntartva.
 Copyright © 2020 Becton, Dickinson and Company. Všechna práva vyhrazena.
 Telif hakki © 2020 Becton, Dickinson and Company. Tüm hakları saklıdır.
 Copyright © 2020 Becton, Dickinson and Company. 保留所有權利。
 저작권 © 2020 Becton, Dickinson and Company. 모든 권리 보유.
 Авторские права © 2020 Компания Becton, Dickinson and Company. Все права защищены.



BD, the BD logo, ClearStream, everlinQ, WavelinQ and all other trademarks are property of Becton, Dickinson and Company or an affiliate
 BD, le logo BD, ClearStream, everlinQ, WavelinQ et toutes les autres marques commerciales appartiennent à Becton, Dickinson and Company ou l'une de ses filiales
 BD, das BD-Logo, ClearStream, everlinQ, WavelinQ und alle anderen Marken sind das Eigentum von Becton, Dickinson and Company oder einer Tochtergesellschaft
 BD, il logo BD, ClearStream, everlinQ, WavelinQ e tutti gli altri marchi commerciali appartengono a Becton, Dickinson and Company o a una società affiliata
 BD, el logotipo de BD, ClearStream, everlinQ, WavelinQ y todas las demás marcas comerciales son propiedad de Becton, Dickinson and Company o de una empresa filial
 BD, het BD-Logo, ClearStream, everlinQ, WavelinQ en alle andere handelsmerken zijn eigendom van Becton, Dickinson and Company of een gelieerde onderneming
 BD, o logótipo BD, ClearStream, everlinQ, WavelinQ e todas as outras marcas comerciais são propriedade da Becton, Dickinson and Company ou de uma afiliada
 Το BD, το λογότυπο BD, ClearStream, τα everlinQ, WavelinQ και όλα τα άλλα εμπορικά σήματα αποτελούν ιδιοκτησία της Becton, Dickinson and Company ή συνδεδεμένων εταιρειών της
 BD, BD-logotet, ClearStream, everlinQ, WavelinQ och alle andre varemærker tilhører Becton, Dickinson and Company eller dets søsterselskaber
 BD, BD-logotypen, ClearStream, everlinQ, WavelinQ och alla övriga varumärken är varumärken och/eller registrerade varumärken som tillhör Becton, Dickinson and Company eller dess dotterbolag
 BD, BD-logo, ClearStream, everlinQ, WavelinQ ja kaikki muut tavaramerkit ovat Becton, Dickinson and Company tai sen tytäryhtiön omaisuutta
 BD, BD logo, ClearStream, everlinQ, WavelinQ og alle andre varemærker tilhører Becton, Dickinson and Company eller et konsernselskap
 BD, logo BD logo, ClearStream, everlinQ, WavelinQ oraz wszystkie inne znaki handlowe są własnością firmy Becton, Dickinson and Company lub jej podmiotów powiązanych
 A BD, a BD logó, ClearStream, az everlinQ, a WavelinQ és minden egyéb védjegy a Becton, Dickinson and Company vagy leányvállalata tulajdona
 BD, logo BD, ClearStream, everlinQ, WavelinQ a veškeré další ochranné známky jsou vlastnictvím společnosti Becton, Dickinson and Company nebo přidružených společností
 BD, BD logosu, ClearStream, everlinQ, WavelinQ ve diğer tüm ticari markalar Becton, Dickinson and Company şirketine veya bağlı kuruluşlarına aittir
 BD, BD 標誌, ClearStream, everlinQ, WavelinQ 與其他所有商標均屬 Becton, Dickinson and Company 或其附屬機構之財產
 BD, BD 로고, ClearStream, everlinQ, WavelinQ 및 기타 모든 상표는 Becton, Dickinson and Company 또는 계열회사
 BD, логотип BD, ClearStream, everlinQ, WavelinQ и все остальные товарные знаки являются собственностью компании Becton, Dickinson and Company или ее аффилированных компаний.

Distributed by:
Bard Peripheral Vascular, Inc.
www.bardpv.com/globalcontact
TEL: 1-800-321-4254

 **Manufacturer:**
ClearStream Technologies, Ltd.
Moyné Upper
Enniscorthy
Co. Wexford, Ireland

Assembled in Ireland

CE
0344