

# An Attractive Non-Surgical AV Fistula Creation Option

The WavelinQ™ EndoAVF System uses two thin, flexible, magnetic catheters and a burst of RF energy to create an AV fistula without open surgery. This non-surgical approach helps to reduce surgical scarring and can minimize arm disfigurement often associated with traditional AV fistula surgery.

**WavelinQ™**  
EndoAVF System



The WavelinQ™ EndoAVF System should not be used in patients who have known central venous stenosis or upper extremity venous occlusion on the same side as the planned AVF creation, who have a known allergy or reaction to any drugs, or who have known adverse reactions to moderate sedation and/or anesthesia.

The WavelinQ™ EndoAVF System is intended for the cutting and coagulation of blood vessel tissue in the peripheral vasculature for the creation of an arteriovenous fistula used for hemodialysis.

Please consult product labels and instructions for use for indications, contraindications, hazards, warnings and precautions.

[crbard.com/peripheral-vascular](http://crbard.com/peripheral-vascular) | [bd.com](http://bd.com)  
[WavelinQSupportEMEA@bd.com](mailto:WavelinQSupportEMEA@bd.com)

BD Switzerland Sarl, Terre Bonne Park – A4, Route De Crassier, 17, 1262 Eysins, Vaud, Switzerland.

© 2021 BD. BD, the BD Logo, and WavelinQ are trademarks of Becton, Dickinson and Company or its affiliates. All Rights Reserved. BD-38933

CE  
0344

