



WavelinQ™

EndoAVF System

Carbon Fiber Arm Board and Fixation Straps

The WavelinQ™ Arm Board is designed to secure and orientate the patient's right or left arm during the WavelinQ™ EndoAVF procedure. It is lightweight, ultra-durable and invisible under fluoroscopy. The arm board and straps help to restrict unwanted movement during the procedure.



Securing the arm for the WavelinQ™ EndoAVF Procedure

Step 1



Step 2



Step 3



Step 4



Step 5



Ordering Information

WavelinQ™ Carbon Fiber Arm Board Item Code

CZ-400-TVA

58.2 cm
(22.9 in)



31.5 cm
(12.4 in)

112.8 cm
(44.4 in)

<8 pounds

Disposable Straps Item Code

TVA-MC-2



Each pack contains three straps: two 2" x 24" and one 2" x 18". Sold as a case of 10 packs.

PHYSICIAN NAME

PHYSICIAN SIGNATURE

REPRESENTATIVE'S NAME

CONTACT PHONE NO.